

# **AccuVax**<sup>®</sup> Vaccine Management System

The First Fully-Automated Vaccine Storage and Handling System that Manages Inventory, Reporting and Ordering for you.

## Efficient and Effective Handling of Vaccine

- Provides assurance that vaccine is at optimal potency when administered
- Improves record keeping accuracy and saves time

## Consistent, Accurate Temperature Management

- 8 data loggers recording temperature and logging 24/7
- Stores vaccine at the exact temperature needed, whether refrigerated or frozen
- Doorless unit maintains stable temperature and protects against temperature excursions, even during power outages

## Built-in Power Backup

- Battery provides at least 9+ hours of coverage in case of power failure – 99% of power outages in the US are shorter than 8 hours
- Provides time to implement emergency plans if necessary

## Stock Rotated Automatically

- Ensures doses nearest to expiry are used first, with no staff time spent on manual stock rotation
- Automatically tracks lot-code and expiration date

## Instant Inventory Reporting

- Full reporting and analytics: inventory, usage patterns, public or private doses
- Helps decision making on inventory levels; tracks public doses separately, generates reports



# AccuVax<sup>®</sup>

Vaccine Management System

The First Fully-Automated Vaccine Storage and Handling System that Manages Inventory, Reporting and Ordering for you.

## Ordering Your Way

- Automated re-ordering (optional), or low inventory prompt
- Customizable ordering compatible with your current process, or fully automated through electronic data transfer

## Proactive Alerting Suite

- Avoid product expiration, stock depletion, temperature excursion
- Instant notification of battery activation
- Alerts sent via text and email

## Remarkably Cost-Effective and Simple to Use

- Intuitive touch screen interface
- No time-consuming training needed – you will be up and running on Day One

## Technical Specifications:

- Capacity: 684 refrigerated doses, 72 frozen
- Size: 28" wide, 28" deep, 56" high
- Communications: LAN, Wi-Fi, and Cellular
- HIPAA compliant
- Integrated 2D barcode reader
- Secure access: fingerprint log-in



For complete information visit [www.AccuVax.com](http://www.AccuVax.com) or call 1-844-TruMed-1